



MODEL 505
PAT 10-2340561

Developed the world's first automatic starter.
Helps to strengthen immunity by generating microcurrent



www.tosung.kr

Made In Korea
TOSUNG CORP

Product specification

- Automatic start-up by installing clutch
 - The three primary colors of the rotor are produced during operation to arouse interest
 - Self-luminous LED lighting provides movement dynamism(Model 707)
 - Eddy current is automatically conducted to help restore bioelectricity
 - Helps to strengthen the immune system of the immunocompromised
 - Model 707 conductive plate 99% pure copper
 - Model 808 conductive plate 99% pure silver
- When using, make sure your hand touches the conductive plate at the bottom of the product.

Size 7.3 x 7.3 x 5.2cm

Product use

- Children's motor development and brain promotion device
- Wrist stretching device (housework, driving, computer)
- Portable wrist exercise equipment
- Mechanism for preventing and improving carpal tunnel syndrome from computer work
- Adult leisure items and strength restoration equipment
- Silver layer cognitive function improvement and dementia prevention device
- Physical therapy equipment for rehabilitation patients who need low-intensity exercise
- Immune strengthening device for immunocompromised persons

OPEN

TOSUNG TOP
Wrist Exercise Goods



MODEL 505
PAT 10-2340561



- ① While holding the front and back sides of the product with your fingers, turn your wrist left and right to rotate the rotor.
- ② Rotate your wrist to the left and right, then bring it back to the original position, and repeat the motion.
- ③ Take turns holding the product with both hands and rotate your wrist repeatedly.
The higher the rotation speed, the brighter the LED and you can feel the dynamics of the movement.
- ④ It is equipped with a magnet, so the moment you hold the product, it generates a microcurrent by eddy current.
This microcurrent can restore the human body's bioelectricity (70~90mV) and increase immunity.
(German Nobel Prize laureate Dr. Erwin Neher won the Nobel Prize for discovering that there are ion channels in human cells, which are regulated by microcurrents)
- ⑤ Models 808 and 909 are made of strong magnets and pure silver conductive material, which can shorten the time to strengthen immunity.
- ⑥ If you are weak, consult your doctor before use.
- ⑦ If this product is arbitrarily disassembled, A/S is not possible.
- ⑧ As magnetic force is generated, avoid contact with precision instruments and store safely.



Wrist Exercise Goods
TOSUNG TOP

www.tosung.kr